










Stapling Product Reference Chart


Endocutters*†

Description	Device	Device Code	Shaft Length	Staple Line Length	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Closed Staple Height	Qty/box
Echelon Flex™ Powered Plus 60mm Stapler		PCEE60A	280mm	60mm	12	3	GST60W	6	<input type="checkbox"/> Vascular/Thin	1.0mm	12
		PSEE60A	340mm				GST60B		<input checked="" type="checkbox"/> Regular	1.5mm	
		PLEE60A	440mm				GST60D		<input checked="" type="checkbox"/> Regular/Thick	1.8mm	
							GST60G		<input checked="" type="checkbox"/> Thick	2.0mm	
							GST60T		<input checked="" type="checkbox"/> Very Thick	2.3mm	
Echelon Flex™ 60mm Stapler		SC60A	280mm	60mm	12		ECR60M	6	<input checked="" type="checkbox"/> Mesentery/Thin	.75mm	
		EC60A	340mm				ECR60B		<input checked="" type="checkbox"/> Regular	1.5mm	
							ECR60D		<input checked="" type="checkbox"/> Regular/Thick	1.8mm	
		LONG60A	440mm				ECR60G		<input checked="" type="checkbox"/> Thick	2.0mm	
Description	Device	Device Code	Shaft Length	Staple Line Length	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Closed Staple Height	Qty/box
Echelon Flex™ Powered 45mm Stapler		PCE45A	280mm	45mm	12		ECR45M	6	<input checked="" type="checkbox"/> Mesentery/Thin	.75mm	
		PSE45A	340mm				ECR45W		<input type="checkbox"/> Vascular/Thin	1.0mm	
		PLE45A	440mm				ECR45B		<input checked="" type="checkbox"/> Regular	1.5mm	
Echelon Flex™ 45mm Stapler		SC45A	280mm	45mm	12	3	ECR45D	6	<input checked="" type="checkbox"/> Regular/Thick	1.8mm	12
		EC45A	340mm				ECR45G		<input checked="" type="checkbox"/> Thick	2.0mm	
		EC45AL	440mm								
Endopath™ ETS Flex 45mm Stapler		ATS45	340mm	45mm	8		TR45W	6	<input type="checkbox"/> Vascular/Thin	1.0mm	
							6R45B	6	<input checked="" type="checkbox"/> Regular	1.5mm	
							TR45G	4	<input checked="" type="checkbox"/> Thick	2.0mm	
Description	Device	Device Code	Shaft Length	Staple Line Length	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Closed Staple Height	Qty/box
Echelon Flex™ Powered Vascular 35mm Stapler		PVE35A	320mm	35mm	12	3	VASECR35	4	<input type="checkbox"/> Vascular/Thin	1.0mm	12


Intraluminal Circular Staplers*

Description	Device	Device Code	Lumen	Head Diameter	Staple Rows	Open Staple Height	Adjustable Closed Staple Height	Qty/box
Curved Intraluminal Stapler		CDH21A	12.4mm	21mm	2	5.5mm	1.0–2.5mm	3
		CDH25A	16.4mm	25mm				
		CDH29A	20.4mm	29mm				
		CDH33A	24.4mm	33mm				
Endoscopic Curved Intraluminal Stapler		ECS21A	12.4mm	21mm	2	5.5mm	1.0–2.5mm	
		ECS25A	16.4mm	25mm				
		ECS29A	20.4mm	29mm				
		ECS33A	24.4mm	33mm				
Straight Intraluminal Stapler		SDH21A	12.4mm	21mm	2	5.5mm	1.0–2.5mm	
		SDH25A	16.4mm	25mm				
		SDH29A	20.4mm	29mm				
		SDH33A	24.4mm	33mm				

Procedure Stapler*




Description	Device	Device Code	Lumen	Head Diameter	Staple Rows	Open Staple Height	Adjustable Closed Staple Height	Qty/box
PPH Procedure Set		PPH03	24.4mm	33mm	2	4mm	.75–1.5mm	3
Includes Proximate® Hemorrhoidal Circular Stapler 33mm, suture threader, circular anal dilator and purse string anoscope								



Curved Cutter Staplers

Description	Device	Device Code	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Open Staple Height	Closed Staple Height	Qty/box
Contour® Curved Cutter Stapler		CS40B	6	3	CR40B CR40G	4	<div>■ Regular</div> <div>■ Thick</div>	3.5mm 4.7mm	1.5mm 2.0mm	3
		CS40G	6							




* Intraluminal Circular Staplers and Procedure Stapler are not reloadable.



Linear Staplers

Description	Device	Device Code	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Adjustable Closed Staple Height	Qty/box
Proximate® Reloadable Linear Stapler 30mm Dial Compression with Heavy Wire option		TLV30	4	3	TRV30	2	<input type="checkbox"/> Vascular/Thin	1.0mm	12
		TL30	4		TR30	2	<input type="checkbox"/> Vascular/Thin <input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick <input checked="" type="checkbox"/> Very Thick	1.0–2.5mm	
		TLH30	4		TRH30	2	<input type="checkbox"/> Vascular/Thin <input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick <input checked="" type="checkbox"/> Very Thick	1.0–2.5mm	
Proximate® Reloadable Linear Stapler 60mm Dial Compression with Heavy Wire option		TL60	4		TR60	2	<input type="checkbox"/> Vascular/Thin <input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick <input checked="" type="checkbox"/> Very Thick	1.0–2.5mm	
		TLH60	4		TRH60	2	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick <input checked="" type="checkbox"/> Very Thick	1.5–2.5mm	
Proximate® Reloadable Linear Stapler 90mm Dial Compression with Heavy Wire option		TL90	4		TR90	2	<input type="checkbox"/> Vascular/Thin <input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick <input checked="" type="checkbox"/> Very Thick	1.0–2.5mm	
		TLH90	4		TRH90	2	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick <input checked="" type="checkbox"/> Very Thick	1.5–2.5mm	




Description	Device	Device Code	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Closed Staple Height	Qty/box
Proximate® Reloadable Linear Stapler 30mm		TX30V	8	3	XR30V	3	<input type="checkbox"/> Vascular/Thin	1.0mm	12
		TX30B	8		XR30B	2	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	
		TX30G			XR30G				
Proximate® Reloadable Linear Stapler 60mm		TX60B	8		XR60B	2	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	
		TX60G			XR60G				


Linear Cutters

Description	Device	Device Code	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Closed Staple Height	Qty/box
Proximate® Linear Cutter 55mm		TVC55	8	3	TVR55	4	<input type="checkbox"/> Vascular/Thin	1.0mm	12
		TLC55	8		TCR55	4	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	
		TCT55	8		TRT55	4			
Proximate® Linear Cutter 75mm		TLC75	8		TCR75	4	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Regular/Thick <input checked="" type="checkbox"/> Thick	1.5mm 1.8mm 2.0mm	
		TCD75	8		TRD75	4			
		TCT75	8		TRT75	4			
Proximate® Linear Cutter 100mm		TLC10	8		TCR10	4	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Thick	1.5mm 2.0mm	
		TCT10	8		TRT10	4			

Description	Device	Device Code	Max Firings	Qty/box	Reload Code	Staple Rows	Tissue	Selectable Closed Staple Height	Qty/box
Ethicon Endo Surgery™ Linear Cutter 55mm		NTLC55	12	3	SR55	6	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Regular/Thick <input checked="" type="checkbox"/> Thick	1.5mm, 1.8mm or 2.0mm	12
Ethicon Endo Surgery™ Linear Cutter 75mm		NTLC75	12		SR75	6	<input checked="" type="checkbox"/> Regular <input checked="" type="checkbox"/> Regular/Thick <input checked="" type="checkbox"/> Thick	1.5mm, 1.8mm or 2.0mm	

Skin Staplers*

Description	Device	Device Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Size	Qty/box
Proximate® Rotating-Head Skin Stapler		PRR35	35	3.9mm	0.53mm	5.7mm	Regular	6
		PRW35	35	3.9mm	0.58mm	6.9mm	Wide	
Proximate® PX Fixed-Head Skin Stapler		PXR35	35	3.9mm	0.53mm	5.7mm	Regular	
		PXW35	35	3.9mm	0.58mm	6.9mm	Wide	
Proximate® Plus MD Fixed-Head Stapler		PMR35	35	3.9mm	0.53mm	5.7mm	Regular	
		PMW35	35	3.9mm	0.58mm	6.9mm	Wide	

Description	Device	Device Code	Qty/box
Squeeze Handle Staple Extrator		PSX	12

Squeeze Handle Staple Extrator is compatible with all Proximate skin staples

* Skin staplers not reloadable